

Fluorescence Measurements: Horiba Jobin Yvon instrument, in Bill Cooper's lab

Calibration.

- Start **FluorEssence** program

- Click on **M**: system will initialize (1 min)

- Once Main experiment menu appears, click on Spectra, then "emission"

- Next, "load": select "water_Raman_FItoolbox.xml" and open.

- click on detector icon (left hand side) then remove S1 and R1 on the formula menu (highlight and then click "remove") Need only S1/R1

Fill cell with distilled water then place in holder (always position cuvette with same side facing front)

Click **Run**

Window will appear at end of scan, asking for the project name. Click **cancel**

Go to File, Export, ASCII data, name the file. Save. Make sure all boxes are unchecked. Click OK.

[Samples can also be calibrated from blank (ie full 2d scan of distilled water) via R-scripts. In this case, however, units will be arbitrary.]

Sample measurement

M

Choose 3D, for 3D scan

-> load method ("3D_EEM.xml")

-> Run

Signal: save S1/R1

project name: cancel

Activate S1/R1 window

-> File -> export -> Ascii data

when exporting, do not include column names etc

once saved, close all plots before taking next sample

This procedure generates 3D excitation-emission plots and corresponding matrices of data. The plots and data are calibrated in "Quinine Sulfate Equivalent" units (QSE). Raman and Raleigh scattering are removed.

Processing via Matlab script "EEM_processing.m"

(change file name at top)

* save in the directory where the data is *

need water raman file

Fluorescence toolbox (Matlab)

Import DataMax File -> Import 53 files *.prn (only need to import first file)

choose ie 5 nm wavelength slit

start: 240 nm excitation

interval: 5 nm

water Raman

-> find file (no extension, edit file name if necessary)

automatically calculates Raman peak]

(leave calibration set)

-> **Proceed**

can then export jpg picture or matrix of data.