## Daily Report: Reaction of Soluble Fe(III) with Catechol (CA) at pH=3

Experiment Date: 11 June, 2019

**Henry Chin** 

## Goal:

- Reproduce the data from Slikboer et al. (2015) with the new chemicals we just purchased
- Serves as a control experiment for follow up work

Molar ratios: Fe: CA [2:1]

#### **Chemicals:**

- KCl Sigma-Aldrich; ACS reagent 99.0-100.5% (74.55 g/mol)
- FeCl<sub>3</sub>.6H<sub>2</sub>O Sigma-Aldrich Iron(III)hexahydrate; ACS reagent 97% (270.30 g/mol)
- C<sub>6</sub>H<sub>6</sub>O<sub>2</sub> Catechol Sigma-Aldrich (110.1 g/mol)
- KOH: BDH 97%; ACS reagent
- HCI: 6.0 N, EMD
- Nanopure water.

## **Procedure:**

#### 1- Solutions preparation:

Solutions of <u>Fe and CA</u> were prepared in 0.01M KCl solution to give the following concentrations after mixing (total volume of reaction mixture=19.5ml): Fe: 2 mM, CA: 1 mM.

#### -Concentrations before mixing:

- 0.01M KCI [v=250 ml, Mwt: 74.55g/mol, mass: 0.1870g]. pH adjusted to 3.01 using HCl.
  - 0.1870/74.55/0.25 = **0.01 M**
- 1.06mM CA [v=100 ml, Mwt: 110.10/mol, mass: 0.0118g]. (using KCl solution)
  - 0.0118/110.10/0.10 = **1.07 mM**
- 80 mM FeCl₃ [v=10 ml, Mwt: 270.03g/mol, mass: 0.2174g]. (using KCl solution)
  - 0.2174/270.03/0.01 = **80.5** mM

#### -Concentrations after mixing:

- (1.06 mM CA × 19ml) ÷ 19.5ml = 1 mM
- $(80 \text{ mM} \text{ Fe} \times 0.5 \text{ml}) \div 19.5 \text{ml} = 2 \text{ mM}$

Rxns

Initial [Fe] = 80 Mm
Initial  $V_{(Fe)}$ = 0.5 ml
After mixing, [Fe] = 2 mM
Initial [CA] = 1.06 Mm
Initial  $V_{(CA)}$ = 19.0 ml
After mixing, [CA] = 1 mM
Molar ratio Fe: CA
2:1

#### 2- Reaction solution:

A) In to 20 mL vial a 19 mL of 1.06 mM catechol was added. The pH was adjusted to 3.03 by HCl solution. Then the vial was wrapped by Al-foil and a magnetic stir bar was added. The reaction started by adding 0.5 mL of 80 mM FeCl<sub>3</sub> to the above CA solutions. Photos were taken after 5 min, 1 hour and 2 hours. UV-Vis spectra were taken after 5 min, 1 hour, and 2 hours. The solution was then syringe-filtered and let dry over night.

## **Data & Results:**


**Table (1):** pH values for prepared solutions **before** starting the reaction (mixing) and their photos:

Solution	KCI	KCl+CA	KCl+FeCl3
Photo	N/A		
рН	3.03	3.11 -> 3.02	1.86

Table (2): pH values and photos for reaction mixtures. Fe: CA [2:1]


Photo			
рН	2.58	2.41	2.38


Table (3): Syringe filter solution and the results:



<sup>\*</sup>Notes: Some Solution was used for UV-Vis and discarded before the process of syringe filtration. This Yield is an underestimation.

## **UV-Vis Spectra:**





(zoom-in)

# **Conclusions and next steps:**

- Experiments successfully reproduced the previously published data
- Next: repeat the same experiment but in the presence of irradiation